

Agilent 6400 Series QQQ LC/MS Techniques and Operation

Course Number R1893A

Volume I

Student Manual





Agilent Technologies

Agilent 6400 Series QQQ LC/MS Techniques and Operation

Course Number R1893A

Volume I

Student Manual

Data Acquisition B.02.01; Qual B.02 SP3; Quant B.03.02

Printed in August, 2009

Notice

The information contained in this document is subject to change without notice.

Agilent Technologies makes no warranty of any kind with regard to this material, including but not limited to the implied warranties of merchantability and fitness for a particular purpose.

Agilent Technologies shall not be liable for errors contained herein or for incidental, or consequential damages in connection with the furnishing, performance, or use of this material.

No part of this document may be photocopied or reproduced, or translated to another program language without the prior written consent of Agilent Technologies, Inc.

Agilent Technologies, Inc
3750 Brookside Pkwy
Suite 100
Alpharetta, GA 30022

© 2009 by Agilent Technologies, Inc.

All rights reserved

Printed in the United States of America

Table of Contents

TRIPLE QUADRUPOLE (QQQ) PRINCIPLES.....	1
In This Section	2
Mass Analysis Basics.....	3
Mass-to-Charge Ratio	4
Performance Characteristics of Mass Analyzers	5
Accuracy	6
Quadrupole Operation.....	7
Quadrupole Mass Filter.....	8
SIM and Scan	11
QQQ vs. other LC/MS Techniques.....	12
Performance Characteristics of Mass Analyzers	13
Tandem Mass Analysis	14
CID.....	16
Tandem Mass Analysis.....	18
Tandem MS in Space	20
MS2 Scan	22
Product Ion Scan	23
Precursor Ion Scan	24
Neutral Loss Scan	25
Single Reaction Monitoring.....	26
Multiple Reaction Monitoring	27
QQQ Applications	29
AGILENT LC/MS QQQ HARDWARE OVERVIEW	31
In This Section	32
QQQ Specifications	33
QQQ System Diagram	34
Ion Sources.....	36
Ion Optics.....	38
Capillary.....	39
Octopole Ion Guide.....	42
MS1	46
Collision Cell	48

MS2.....	53
HED	54
Vacuum System	57
MassHunter WorkStation.....	63
Compliance Software - Optional	65
IONIZATION SOURCES AND JET STREAM TECHNOLOGY.....	67
In This Section	68
QQQ Ion Sources.....	69
LC/MS Techniques	70
Electrospray	71
Electrospray Ion Suppression	82
ESI with Jet Stream Technology.....	89
APCI	96
Multimode.....	103
HPLC-Chip Cube.....	111
HPLC-Chip MS – Acquisition Software	114
APPI.....	115
Laboratory.....	120
AGILENT TRIPLE QUADRUPOLE TUNING.....	121
In This Section	122
QQQ Tuning	123
QQQ Checktune.....	136
QQQ Autotune	143
QQQ Tuning	159
Tune Troubleshooting.....	160
Multimode Source QQQ Tuning	164
Manual Tuning.....	165
Agilent 6460 Tuning.....	166
LABORATORY EXERCISE: AGILENT QQQ ROUTINE MAINTENANCE AND TUNING.....	181
Preparation:.....	182
Checktune and Autotune.....	183
Routine Maintenance Checklist.....	184
IMPORTANT HPLC TIPS AND MAINTENANCE	189

In This Section	190
Fittings	191
Capillaries	193
Extra-Column Volume.....	194
General Agilent 1200 Recommendations	195
Solvent Preparations	198
Cleaning the HPLC.....	202
Solvent Miscibility.....	204
Solvent Preparations	205
Keys to Success.....	207
Solvent Compressibility.....	208
Preparing the HPLC	210
Sample Preparation	211
HPLC Columns.....	217
Agilent 1200 HPLC Maintenance Areas	220
Recommended Tests	222
Pump Maintenance.....	223
Is My Pump Well?	226
Autosampler Maintenance	227
Thermostatted Column Compartment.....	230
Pressure Problems.....	234
METHODS AND COLUMNS FOR LC/MS.....	239
In This Section	240
Sample Information and Considerations.....	241
Selecting Columns and Conditions for LC/MS	242
Adapting Existing LC Methods	244
Compatible Solvents	245
Mobile Phase Buffers.....	246
Bonded Phase Selection.....	247
Column Configuration	252
Particle Size	253
Rapid Resolution Columns	258
Available Columns.....	260

Comparison Standard and Rapid Resolution Columns.....	261
Increasing Temperature to Reduce Analysis Time.....	262
Rapid Resolution Column Hardware.....	263
Method Transfer to Rapid Resolution Columns.....	264
Superficially Porous Particles.....	266
Specialized Instrumentation.....	270
COLLISION INDUCED DISSOCIATION (CID).....	277
In This Section.....	278
CID.....	279
Interpretation of CID Spectra.....	281
Site of Ionization.....	283
Common Neutral Loss.....	284
Ion Stability.....	285
Fragmentation Mechanisms.....	286
Single Bond Cleavage.....	287
Cyclizations.....	288
Multiple Cleavages.....	289
Ring Opening.....	290
Peptide Sequencing.....	291
Charge Remote Fragmentation.....	292
Fragmentation Example.....	294
QQQ METHOD DEVELOPMENT AND OPTIMIZATION.....	299
In This Section.....	300
Method Development Overview.....	301
Source Parameters.....	303
Agilent 6460 Parameters.....	304
Confirming Precursor Ion.....	309
Parameter Optimization.....	310
Flow Injection Analysis.....	313
MassHunter Optimizer.....	321
Infusion.....	338
Summary.....	339
QQQ DATA ACQUISITION.....	341

In This Section	342
Data Acquisition	343
Instrument Status	345
Instrument Actuals	346
Instrument Control	350
Standby	351
Single Sample Method.....	352
Method Pane	353
Autosampler Parameters	354
Injector Program	355
Binary Pump Parameters.....	356
Column Compartment.....	359
DAD Parameters	360
QQQ Parameters	362
Agilent 6460 Parameters.....	380
QQQ Parameters	381
Single Injection	387
Real-time Plot	388
Connections with Data Analysis.....	390
LABORATORY EXERCISE: METHOD DEVELOPMENT AND OPTIMIZATION	393
Preparation	394
Prepare the QQQ and Agilent 1200 LC.....	395
Conduct an MS2 Scan Analysis.....	396
Use FIA to Optimize Fragmentor Voltage and Collision Energy.....	399
Using MassHunter Optimizer to Determine Optimum Fragmentor Voltage and Collision Energy.....	403

Agilent 6400 Series QQQ LC/MS Techniques and Operation

Course Number R1893A

Volume II

Student Manual





Agilent Technologies

Agilent 6400 Series QQQ LC/MS Techniques and Operation

Course Number R1893A

Volume II

Student Manual

Data Acquisition B.02.01; Qual B..02 SP3; Quant B.03.02
Printed in September, 2009

Notice

The information contained in this document is subject to change without notice.

Agilent Technologies makes no warranty of any kind with regard to this material, including but not limited to the implied warranties of merchantability and fitness for a particular purpose.

Agilent Technologies shall not be liable for errors contained herein or for incidental, or consequential damages in connection with the furnishing, performance, or use of this material.

No part of this document may be photocopied or reproduced, or translated to another program language without the prior written consent of Agilent Technologies, Inc.

Agilent Technologies, Inc
3750 Brookside Pkwy
Suite 100
Alpharetta, GA 30022

© 2009 by Agilent Technologies, Inc.

All rights reserved

Printed in the United States of America

Table of Contents

WORKLISTS FOR AUTOMATION	1
In This Section	2
Worklist Pane.....	3
Create Acquisition Methods	4
New Worklist or Worklist from Template.....	5
Worklist Columns	6
Add Samples	8
Start Worklist	14
Partial Worklist	16
Ways to Create Worklists	17
Worklist Run Parameters	24
Scripts	27
Worklist Reports	32
Save Worklist.....	35
LABORATORY EXERCISE: WORKLISTS	37
Preparation	38
Preparing the LC/MSD QQQ.....	39
Worklists	40
BASIC QUALITATIVE DATA ANALYSIS	43
In This Section	44
Open Qualitative Data Analysis.....	45
Configuration	46
Data Analysis Layout.....	47
Compound Labels	49
Method Explorer	51
Methods.....	55
Open Data Files.....	57
Open Data Files.....	59
Opening Multi-Signal Files	60
Zooming Data	61
Chromatogram Functions.....	62

Chromatogram Functions.....	63
Anchoring	68
Integration Functions	69
System Suitability	75
Spectrum Functions	76
Spectrum Display Options	77
Spectrum Functions	78
Reports	92
Find Compounds by Formula	94
Tools	99
LABORATORY EXERCISE: BASIC QUALITATIVE DATA ANALYSIS	103
In This Exercise, You Will Need:.....	104
Open the Qualitative Application and Restore Defaults.....	105
Opening Single or Multiple Data Files	106
Zooming In and Zooming Out	107
Displaying Spectra	108
Chromatogram Extraction.....	112
Signal-to-Noise Calculation.....	114
Find Compounds by Formula	116
Exploring Different Workflows.....	117
Close Data Files	118
QUALITATIVE ANALYSIS AUTOMATION	119
In This Section	120
Data Analysis Automation.....	121
File Open Actions	122
Worklist Actions	131
Compound Automation Steps	132
Automated Data Analysis	133
LABORATORY EXERCISE: QUALITATIVE ANALYSIS AUTOMATION	135
Preparations:	136
File Open Actions	137
Set Actions to Run When Opening a Data File	139

ASSEMBLING THE MRM METHOD.....	141
In This Section	142
Overview.....	143
MS2 Scan Method.....	144
MS2 Scan to MRM Method.....	145
Import Compound Data from Optimizer	147
Dwell Time	151
Calibrating the MRM Method	153
View in Qualitative Analysis.....	154
Use Walk the Chromatogram.....	155
Time Segments.....	156
Dynamic MRM Scan Type	157
Delta EMV	163
Optimizing Chromatography	164
Calibrating the MRM Method	165
Fine Tuneing the MRM Method.....	166
Evaluating the MRM Method.....	167
Method Validation	168
LABORATORY EXERCISE: ASSEMBLING THE MRM METHOD.....	171
Preparation	172
Assemble the MRM Method.....	173
QQQ QUANTIFICATION: BUILD A QUANTIFICATION METHOD	177
In This Section	178
Basic Workflow	179
Open MassHunter Quantification	180
Batch	181
Add Samples to Batch.....	182
Develop Method using Calibration Standard.....	183
Open Method Editing.....	184
New Method Easy Development	185
MRM Compound Setup: Information Added Automatically from Acquisition Method .	186
New Method Manual Setup	187
Method Edit Tasks	188

Retention Time Setup	189
ISTD Setup: 4 columns.....	194
Concentration Setup.....	195
Create Levels	196
Copy Calibration Levels	197
Setup Differing Concentrations	198
Qualifiers.....	199
Calibration Curve.....	200
Global Setup.....	201
Validate Quantification Method	202
Save and Exit	203
Return to Quantitative Analysis Main Screen	204
Analyze Batch.....	205
After Analyze Batch is Executed.....	206
Compound Information Display Options.....	207
Messages and Outliers	208
Calibration Curve.....	210
Remove Calibration Data Points.....	211
Curve Fit Assistant.....	212
Sort and See Confidence Interval	217
Save Batch with Results.....	218
LABORATORY EXERCISE: BASIC QUANTITATIVE ANALYSIS – TRIPLE QUADRUPOLE.....	219
Preparations.....	220
Create the Batch.....	221
Create the Quantitative Method	223
Analyze the Batch	229
Adjust the Curve Fit.....	231
Navigate Batch-at-a-Glance.....	233
ADVANCED QQQ QUANTIFICATION	245
In This Section	246
Batch Table – Key Information	247
Accuracy Column with Outlier Flags	248
When the Accuracy Column Has Outliers.....	249

When the Accuracy Column has Outliers.....	250
Zoom In and View Calibration Curve	253
Full View of Calibration Curve	254
After Re-Analyzing the Batch	255
Next Problem: The Blank	256
Calibration Curve.....	257
Calibration Curve: Significant Positive x Intercept.....	258
Calibration Curve: Significant + x Intercept.....	261
Other Important Columns in the Batch Table: Qualifier Results	267
Procedure for Updating Qualifier Ion Ratios.....	268
Concentrations in Batch Table.....	269
Adding Batch Table Columns.....	270
Adding Batch Columns.....	272
Other Important Columns in the Batch Table: ISTD Response	275
Advanced QQQ Quantification: Summary.....	276
LABORATORY EXERCISE: MORE QUANTITATION.....	279
Preparations.....	280
Integrate Without Parameters	281
Display Outliers	288
LABORATORY EXERCISE: ADVANCED QUANTITATIVE ANALYSIS II.....	295
Preparations.....	296
Create a New Batch and Refine the Calibrations.....	297
Evaluating the Calibration Curve.....	304
MASSHUNTER REPORTING.....	313
In This Section	314
MassHunter Reporting Fundamentals.....	315
Printing Reports	318
Exporting Reports	322
Report Generation Summary	323
Customized Report Templates.....	324
LABORATORY EXERCISE: MASSHUNTER REPORTING.....	335
Preparations.....	336
Generate Reports.....	337

Customizing the Report Template	339
QQQ MAINTENANCE AND TROUBLESHOOTING.....	351
In This Section	352
Reduce Instrument Downtime	353
Common Causes Instrument Failure.....	354
Typical Background.....	355
General Tips.....	356
Objectives	376
Preventive Maintenance Checklist.....	377
Troubleshooting	379